


Product Datasheet

NONO antibody (orb1321849)

Description	Mouse monoclonal antibody to NONO
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 184-385 of human NONO(NP_031389) produced in E.coli.
Target	NONO
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	54.1 kDa
Uniprot ID	Q15233
NCBI	NM_007363
Dilution Range	WB 1:200~1000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...