


Product Datasheet

NQO2 antibody (orb1321866)

Description	Mouse monoclonal antibody to NQO2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human NQO2(NP_000895) produced in HEK293T cell.
Target	NQO2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	25.7 kDa
Uniprot ID	P16083
NCBI	NM_000904
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...