


Product Datasheet

ACSF2 antibody (orb1321898)

Description	Mouse monoclonal antibody to ACSF2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human ACSF2(NP_079425) produced in HEK293T cell.
Target	ACSF2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	67.9 kDa
Uniprot ID	Q96CM8
NCBI	NM_025149
Dilution Range	IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...