


Product Datasheet

CRYZL1 antibody (orb1321932)

Description	Mouse monoclonal antibody to CRYZL1
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human CRYZL1(NP_665857) produced in HEK293T cell.
Target	CRYZL1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	38.5 kDa
Uniprot ID	O95825
NCBI	NM_145858
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...