


Product Datasheet

ATRIP antibody (orb1321945)

Description	Mouse monoclonal antibody to ATRIP
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 42-309 of human ATRIP(NP_569055) produced in E.coli.
Target	ATRIP
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	85.7 kDa
Uniprot ID	Q8WXE1
NCBI	NM_130384
Dilution Range	WB 1:500~2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...