


Product Datasheet

C6orf64 (SAYSD1) antibody (orb1321952)

Description	Mouse monoclonal antibody to C6orf64 (SAYSD1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human C6orf64(NP_060792) produced in HEK293T cell.
Target	SAYSD1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	20 kDa
Uniprot ID	Q9NPB0
NCBI	NM_018322
Dilution Range	WB 1:500~2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...