


Product Datasheet

Ribophorin II (RPN2) antibody (orb1321976)

Description	Mouse monoclonal antibody to Ribophorin II (RPN2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human RPN2(NP_002942) peoduced in HEK293T cell.
Target	RPN2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	67 kDa
Uniprot ID	P04844
NCBI	NM_002951
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...