


Product Datasheet

FXYD3 antibody (orb1321984)

Description	Mouse monoclonal antibody to FXYD3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human FXYD3(NP_005962) produced in HEK293T cell.
Target	FXYD3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	7.1 kDa
Uniprot ID	Q14802
NCBI	NM_005971
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...