

Product Datasheet

Transferrin Receptor 2 (TFR2) antibody (orb1321989)

Description Mouse monoclonal antibody to Transferrin Receptor 2 (TFR2)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human TFR2(NP_003218) produced in HEK293T cell.

Target TFR2

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 88.6 kDa

Uniprot ID [Q9UP52](#)

NCBI [NM_003227](#)

Dilution Range WB 1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...