


Product Datasheet

TAF2 antibody (orb1321991)

Description	Mouse monoclonal antibody to TAF2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 782-1121 of human TAF2(NP_003175) produced in E.coli.
Target	TAF2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	136.8 kDa
Uniprot ID	Q6P1X5
NCBI	NM_003184
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...