


Product Datasheet

TMEFF2 antibody (orb1322017)

Description	Mouse monoclonal antibody to TMEFF2
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human TMEFF2(NP_057276) produced in HEK293T cell.
Target	TMEFF2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	41.2 kDa
Uniprot ID	Q9UIK5
NCBI	NM_016192
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...