


Product Datasheet

LIM1 (LHX1) antibody (orb1322038)

Description	Mouse monoclonal antibody to LIM1 (LHX1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 100-362 of human LHX1(NP_005559) produced in HEK293T cell.
Target	LHX1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	44.6 kDa
Uniprot ID	P48742
NCBI	NM_005568
Dilution Range	WB 1:500~1000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...