

Product Datasheet

Phosphoribosyl pyrophosphate amidotransferase (PPAT) antibody (orb1322048)

Description Mouse monoclonal antibody to Phosphoribosyl pyrophosphate amidotransferase (PPAT)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Human recombinant protein fragment corresponding to amino acids 42-278 of human PPAT(NP_002694) produced in E.coli.

Target PPAT

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 57.2 kDa

Uniprot ID [Q06203](#)

NCBI [NM_002703](#)

Dilution Range WB 1:200 - 1:1000

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...