


Product Datasheet

CD10 (MME) antibody (orb1322049)

Description	Mouse monoclonal antibody to CD10 (MME)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human MME(NP_009218) produced in HEK293T cell.
Target	MME
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	85.3 kDa
Uniprot ID	P08473
NCBI	NM_007287
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...