


Product Datasheet

SCHIP1 antibody (orb1322061)

Description	Mouse monoclonal antibody to SCHIP1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human SCHIP1(NP_055390) produced in HEK293T cell.
Target	SCHIP1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	53.3 kDa
Uniprot ID	P0DPB3
NCBI	NM_014575
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...