


Product Datasheet

CNOT4 antibody (orb1322091)

Description	Mouse monoclonal antibody to CNOT4
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Human recombinant protein fragment corresponding to ammino acids 190-455 of human CNOT4(NP_037448) produced in E.coli.
Target	CNOT4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	70.2 kDa
Uniprot ID	O95628
NCBI	NM_013316
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...