


Product Datasheet

PTOP (ACD) antibody (orb1322136)

Description	Mouse monoclonal antibody to PTOP (ACD)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human ACD(NP_001075955) produced in HEK293T cell.
Target	ACD
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	57.6 kDa
Uniprot ID	Q96AP0
NCBI	NM_001082486
Dilution Range	WB 1:2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...