


Product Datasheet

Cytochrome P450 2B6 (CYP2B6) antibody (orb1322178)

Description	Mouse monoclonal antibody to Cytochrome P450 2B6 (CYP2B6)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human CYP2B6(NP_000758) produced in HEK293T cell.
Target	CYP2B6
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	56.1 kDa
Uniprot ID	P20813
NCBI	NM_000767
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...