


Product Datasheet

TIMP2 antibody (orb1322190)

Description	Mouse monoclonal antibody to TIMP2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TIMP2(NP_003246) produced in HEK293T cell.
Target	TIMP2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21.7 kDa
Uniprot ID	P16035
NCBI	NM_003255
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...