

Product Datasheet

Epsin 2 (EPN2) antibody (orb1322199)

Description	Mouse monoclonal antibody to Epsin 2 (EPN2)
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human EPN2(NP_055779) produced in HEK293T cell.
Target	EPN2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	68.3 kDa
Uniprot ID	O95208
NCBI	NM_014964
Dilution Range	WB 1:500~2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...