


Product Datasheet

TLE2 antibody (orb1322248)

Description	Mouse monoclonal antibody to TLE2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human TLE2(NP_003251) produced in HEK293T cell.
Target	TLE2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	79.7 kDa
Uniprot ID	Q04725
NCBI	NM_003260
Dilution Range	WB 1:2000, FLOW 1:100, IF: 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...