

Product Datasheet

AKT1 antibody (orb1322250)

Description	Mouse monoclonal antibody to AKT1
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human AKT1(NP_005154) produced in HEK293T cell.
Target	AKT1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	55.5 kDa
Uniprot ID	P31749
NCBI	NM_005163
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...