


Product Datasheet

NMT2 antibody (orb1322277)

Description	Mouse monoclonal antibody to NMT2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human NMT2(NP_004799) produced in HEK293T cell.
Target	NMT2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	56.8 kDa
Uniprot ID	O60551
NCBI	NM_004808
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...