

Product Datasheet

GIRK1 (KCNJ3) antibody (orb1322280)

Description	Mouse monoclonal antibody to GIRK1 (KCNJ3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 279-501 of human KCNJ3(NP_002230) produced in E.coli.
Target	KCNJ3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	56.4 kDa
Uniprot ID	P48549
NCBI	NM_002239
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.

