



Product Datasheet

14-3-3 antibody (orb13223)


Description	14-3-3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human 14-3-3 (174-245/245 aa)
Target	14-3-3 Alpha + Beta + Gamma + Delta + Epsilon
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27 kDa
Uniprot ID	P31946
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IF=1:100-500, Flow-Cyt=1µg /test., IHC-F=1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis
of rat brain lysat...



rat heart
lysates(Lane2) using
14-3-3 an...



Immunohistochemical
analysis of formalin...