


Product Datasheet

TPSG1 antibody (orb1322304)

Description	Mouse monoclonal antibody to TPSG1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 20-283 of human TPSG1(NP_036599) produced in E.coli.
Target	TPSG1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	32 kDa
Uniprot ID	Q9NRR2
NCBI	NM_012467
Dilution Range	WB 1:2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...