

Product Datasheet

Cystatin S (CST4) antibody (orb1322335)

Description Mouse monoclonal antibody to Cystatin S (CST4)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human CST4(NP_001890) produced in HEK293T cell.

Target CST4

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 14.1 kDa

Uniprot ID [P01036](#)

NCBI [NM_001899](#)

Dilution Range WB 1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...