


Product Datasheet

PFKP antibody (orb1322432)

Description	Mouse monoclonal antibody to PFKP
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human PFKP(NP_002618) produced in HEK293T cell.
Target	PFKP
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	85.4 kDa
Uniprot ID	Q01813
NCBI	NM_002627
Dilution Range	WB 1:500~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...