


Product Datasheet

PDSS2 antibody (orb1322453)

Description	Mouse monoclonal antibody to PDSS2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human PDSS2(NP_065114) produced in HEK293T cell.
Target	PDSS2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	43.9 kDa
Uniprot ID	Q86YH6
NCBI	NM_020381
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...