


Product Datasheet

DPP9 antibody (orb1322468)

Description	Mouse monoclonal antibody to DPP9
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human DPP9(NP_631898) produced in HEK293T cell.
Target	DPP9
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	96.4 kDa
Uniprot ID	Q86TI2
NCBI	NM_139159
Dilution Range	WB 1:500~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...