


Product Datasheet

TAL1 antibody (orb1322525)

Description	Mouse monoclonal antibody to TAL1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human TAL1(NP_003180) produced in HEK293T cell.
Target	TAL1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	34.1 kDa
Uniprot ID	P17542
NCBI	NM_003189
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...