


Product Datasheet

CENPH antibody (orb1322538)

Description	Mouse monoclonal antibody to CENPH
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human CENPH(NP_075060) produced in HEK293T cell.
Target	CENPH
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	28.3 kDa
Uniprot ID	Q9H3R5
NCBI	NM_022909
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...