


Product Datasheet

NUDT18 antibody (orb1322568)

Description	Mouse monoclonal antibody to NUDT18
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human NUDT18(NP_079091) produced in HEK 293T cell.
Target	NUDT18
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.3 kDa
Uniprot ID	Q6ZVK8
NCBI	NM_024815
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...