

## Product Datasheet

### Adiponectin (ADIPOQ) antibody (orb1322572)

|                            |                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to Adiponectin (ADIPOQ)                                                                           |
| <b>Species/Host</b>        | Mouse                                                                                                                       |
| <b>Reactivity</b>          | Human                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                |
| <b>Tested Applications</b> | FC, IF, IHC, WB                                                                                                             |
| <b>Immunogen</b>           | Human recombinant protein fragment corresponding to amino acids 19-244 of human ADIPOQ(NP_004788) produced in HEK293T cell. |
| <b>Target</b>              | ADIPOQ                                                                                                                      |
| <b>Note</b>                | For research use only                                                                                                       |
| <b>Isotype</b>             | IgG1                                                                                                                        |
| <b>Clonality</b>           | Monoclonal                                                                                                                  |
| <b>MW</b>                  | 24.9 kDa                                                                                                                    |
| <b>Uniprot ID</b>          | <a href="#">Q15848</a>                                                                                                      |
| <b>NCBI</b>                | <a href="#">NM_004797</a>                                                                                                   |
| <b>Dilution Range</b>      | WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                             |

