


Product Datasheet

FAM119A (METTL21A) antibody (orb1322590)

Description	Mouse monoclonal antibody to FAM119A (METTL21A)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human FAM119A(NP_660323) produced in HEK293T cell.
Target	METTL21A
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	24.4 kDa
Uniprot ID	Q8WXB1
NCBI	NM_145280
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...