

Product Datasheet

Adiponectin (ADIPOQ) antibody (orb1322591)

Description	Mouse monoclonal antibody to Adiponectin (ADIPOQ)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 19-244 of human ADIPOQ(NP_004788) produced in HEK293T cell.
Target	ADIPOQ
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	24.9 kDa
Uniprot ID	Q15848
NCBI	NM_004797
Dilution Range	WB 1:500, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.

