

Product Datasheet

Nucleotide binding protein like (NUBPL) antibody (orb1322628)

Description Mouse monoclonal antibody to Nucleotide binding protein like (NUBPL)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, WB

Immunogen Human recombinant protein fragment corresponding to amino acids 1-250 of human NUBPL(NP_079428) produced in E.coli.

Target NUBPL

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 33.9 kDa

Uniprot ID [Q8TB37](#)

NCBI [NM_025152](#)

Dilution Range WB 1:2000, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...