


Product Datasheet

TBC1D21 antibody (orb1322649)

Description	Mouse monoclonal antibody to TBC1D21
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TBC1D21(NP_699187) produced in HEK293T cell.
Target	TBC1D21
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	39 kDa
Uniprot ID	Q8IYX1
NCBI	NM_153356
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...