


Product Datasheet

ZNF397 antibody (orb1322666)

Description	Mouse monoclonal antibody to ZNF397
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human ZNF397(NP_115723) produced in HEK293T cell.
Target	ZNF397
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	30.9 kDa
Uniprot ID	Q8NF99
NCBI	NM_032347
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...