


Product Datasheet

TP73-AS1 antibody (orb1322671)

Description	Mouse monoclonal antibody to TP73-AS1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human KIAA0495(NP_997189) produced in HEK293T cell.
Target	KIAA0495
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	20.4 kDa
NCBI	NM_207306
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...