


Product Datasheet

SIRT6 antibody (orb1322681)

Description	Mouse monoclonal antibody to SIRT6
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 88-352 of human SIRT6(NP_057623) produced in E.coli.
Target	SIRT6
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	38.9 kDa
Uniprot ID	Q8N6T7
NCBI	NM_016539
Dilution Range	WB 1:500~2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...