


Product Datasheet

Ribophorin I (RPN1) antibody (orb1322697)

Description	Mouse monoclonal antibody to Ribophorin I (RPN1)
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human RPN1(NP_002941) produced in HEK293T cell.
Target	RPN1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	66.2 kDa
Uniprot ID	P04843
NCBI	NM_002950
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...