


Product Datasheet

EIF2B3 antibody (orb1322712)

Description	Mouse monoclonal antibody to EIF2B3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human EIF2B3(NP_065098) produced in HEK293T cell.
Target	EIF2B3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	50.1 kDa
Uniprot ID	Q9NR50
NCBI	NM_020365
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...