


Product Datasheet

PIG3 (TP53I3) antibody (orb1322738)

Description	Mouse monoclonal antibody to PIG3 (TP53I3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TP53I3(NP_671713) produced in HEK293T cell.
Target	TP53I3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.4 kDa
Uniprot ID	Q53FA7
NCBI	NM_147184
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...