


Product Datasheet

PRMT2 antibody (orb1322742)

Description	Mouse monoclonal antibody to PRMT2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human PRMT2(NP_001526) produced in HEK293T cell.
Target	PRMT2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	48.9 kDa
Uniprot ID	P55345
NCBI	NM_001535
Dilution Range	WB 1:500, IF: 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...