


Product Datasheet

MTFMT antibody (orb1322761)

Description	Mouse monoclonal antibody to MTFMT
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 175-389 of human MTFMT(NP_640335) produced in E.coli.
Target	MTFMT
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	43.7 kDa
Uniprot ID	Q96DP5
NCBI	NM_139242
Dilution Range	WB 1:500~2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...