


Product Datasheet

Proteasome subunit beta type 4 (PSMB4) antibody (orb1322769)

Description	Mouse monoclonal antibody to Proteasome subunit beta type 4 (PSMB4)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human PSMB4(NP_002787) produced in HEK293T cell.
Target	PSMB4
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	29 kDa
Uniprot ID	P28070
NCBI	NM_002796
Dilution Range	WB 1:500, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...