


Product Datasheet

TAZ (WWTR1) antibody (orb1322784)

Description	Mouse monoclonal antibody to TAZ (WWTR1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human WWTR1(NP_057287) produced in HEK293T cell.
Target	WWTR1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	43.9 kDa
Uniprot ID	Q9GZV5
NCBI	NM_015472
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...