


Product Datasheet

PDE4 (PDE4B) antibody (orb1322788)

Description	Mouse monoclonal antibody to PDE4 (PDE4B)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human PDE4B(NP_002591) produced in HEK293T cell.
Target	PDE4B
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	83.2 kDa
Uniprot ID	Q07343
NCBI	NM_002600
Dilution Range	WB 1:500~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...