


Product Datasheet

CTBP1-DT antibody (orb1322799)

Description	Mouse monoclonal antibody to CTBP1-DT
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-186 of human C4orf42(NP_443093) produced in E.coli.
Target	C4orf42
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	19.5 kDa
NCBI	NM_052861
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...